

# ISA

From Experience comes Excellence



AFI  
Compact centrifuge



## Rotors

Tube	Capacity	Speed
5 mL	40	4700Rpm 3729 x g
15 mL	12	6700 Rpm 6324 x g
50 mL	6	6700 Rpm 6073 x g
100mL	6	6700 Rpm 5320 x g
Microplate	8	4500Rpm 2400 x g

## Caracteristics

### Dimensions (HxWxD)

Ventilated : 34x 39x 51 cm

Refrigerated: 34x 46x 72 cm

### Speed max:

Swing out: 4500 Rpm/ 3396 x g

Angular rotor : 6700Rpm /6324 x g

### AFI Unique points:

- Blue electronic plug & play module
- Digital LCD display, angled for better for visibility
- Touch Start/Stop button accessible with gloved hand
- Interactive communication with the user
- Front USB connection



• Should it be necessary or useful without undermining any model essential features, technical and stylistic characteristics are subject to be changed without manufacturer's previous notice.

M A D E I N F R A N C E



**AFI** Centrifuge

3 Rue Nicolas Copernic - 53200 Château Gontier - France

Tel. +33 2 43 06 15 19

<http://www.afigroups.com>

e-mail: [contact@afigroups.com](mailto:contact@afigroups.com)

*Scan me to discover  
AFI catalogs*

