

LINA

From Experience comes Excellence



AFI centrifuge
Micro & Multifonction



Rotors

Tube	Capacity	Speed
1,5/2mL	12	23 000 Rpm 39 625 x g
	48	14 000 Rpm 21 255 x g
5 mL 90°	8	10 000 Rpm 11 739 x g
10 mL	12	6000 Rpm 4629x g
15mL	12	6000 Rpm 4629x g
50mL	6	8330 Rpm 8766 x g
Microplate	6	4700 Rpm
Deepwell	2	2223 x g

Caracteristics

Dimensions (HxWxD)

Refrigerated: 28 x 35 x 64 cm

Ventilated: 28 x 31 x 40 cm

Speed Max :

Angular rotor: 23 000 Rpm/ 39 625 x g

AFI Unique points:

- Blue electronic plug & play module
- Ultra quiet centrifuge
- Interactive communication
- Illuminated display visible at 7 meters
- Gesture control: start or stop a cycle, open the centrifuge lid without any contact for optimal hygiene.



• Should it be necessary or useful without undermining any model essential features, technical and stylistic characteristics are subject to be changed without manufacturer's previous notice.

M A D E I N F R A N C E



AFI Centrifuge

3 Rue Nicolas Copernic - 53200 Château Gontier - France

Tel. +33 2 43 06 15 19

<http://www.afigroups.com>

e-mail: contact@afigroups.com

Scan me to discover

AFI catalogs

