

# LISA

*From Experience comes Excellence*



**Routine MultiLab**  
**AFI Centrifuge**



## Rotors

Tube	Capacity	Speed
<b>5 mL</b>	<b>92</b>	5100Rpm 5816 <sub>x g</sub>
10mL	76	5100Rpm 5816 <sub>x g</sub>
15mL	60	5100Rpm 5816 <sub>x g</sub>
50 mL	20	5100Rpm 5816 <sub>x g</sub>
<b>100mL</b>	<b>8</b>	<b>14 000 Rpm</b> <b>25 200 x g</b>
<b>625 mL</b>	<b>4</b>	5100Rpm 5816 <sub>x g</sub>
<b>Microplate</b>	<b>16</b>	5100Rpm 5816 <sub>x g</sub>

## Characteristics

### Dimensions (HxLxD)

Ventilated : 34x46x56 cm

Refrigerated: 32x 68x 58 cm

### Max speed:

Swing out rotor: 5100 Rpm/ 5816 x g

Angular rotor: 18 000 Rpm/31 514 xg

### AFI Unique points:

- Blue electronic plug & play module
- Interactive communication with the user
- Click Spin: Tool-free accessory interchangeability system
- Large colour touch screen, gloved hand compatible
- Silent centrifuge



• Should it be necessary or useful without undermining any model essential features, technical and stylistic characteristics are subject to be changed without manufacturer's previous notice.

M A D E I N F R A N C E



**AFI** Centrifuge

3 Rue Nicolas Copernic - 53200 Château Gontier - France  
Tel. +33 2 43 06 15 19  
<http://www.afigroups.com>  
e-mail: [contact@afigroups.com](mailto:contact@afigroups.com)

*Scan me to discover  
AFI catalogs*

