

LOREENA

From Experience comes Excellence



**AFI Routine
MultiLab Centrifuge**

Rotors

Tube	Capacity	Speed
5 mL	200	4250 Rpm 4241 x g
50mL	40	4250 Rpm 4241 x g
100mL	8	13 500 Rpm 23 432 x g
250 mL	6	9000Rpm 13 131 x g
1000 mL	4	4250 Rpm 4241 x g
Microplate	24	4250 Rpm 4241 x g

NEW Microplates adapter

NEW Capacity : 4*1L

Caracteristics

Dimensions (HxWxD):

Bench Refrigerated 35x77x64 cm

Floor model : 83x 57x 68 cm

Speed max :

Swing out: 4600 Rpm/ 4731 x g

Angular rotor: 16 000 Rpm/ 28 198 x g

AFI Unique Points:

- Blue electronic plug & play module
- Interactive communication with the user
- Click Spin: Tool-free accessory interchangeability system
- Large colour touch screen, gloved hand compatible
- Silent centrifuge



• Should it be necessary or useful without undermining any model essential features, technical and stylistic characteristics are subject to be changed without manufacturer's previous notice.

M A D E I N F R A N C E



AFI Centrifuge

3 Rue Nicolas Copernic - 53200 Château Gontier - France

Tel. +33 2 43 06 15 19

<http://www.afigroups.com>

e-mail: contact@afigroups.com

*Scan me to discover
AFI catalogs*

