

# 10 µl Graduated TipOne® Tip, Racked (non-sterile)

TipOne® is the original refillable, reusable and recyclable pipette tip system. Recent improvements to the design of the system, based on customer-feedback, is why TipOne® is one of the leading tips for research labs.

Convenience and quality are at the heart of the system and TipOne® remains the same tip that you have used and trusted over the years.

## Tip fit

One tip fits all! The tip collar is engineered to fit all popular pipettes.

## Manufacture

The use of multiple cavity tools for the injection moulding process guarantees high conformity within one batch. The automated moulding and assembling procedures are carried out under clean-room conditions to eliminate contamination.

## Unique collar

Our newest tips can be recognised by the jagged design on the tip collar.

## Material

All TipOne® tips are manufactured from pure, virgin polypropylene to ensure the highest quality product.

## Green tip design

A reduction of up to 40 % in the material used for the production of TipOne® tips helps to save valuable resources and is beneficial for the environment.

## Highly polished moulds

For an extremely smooth surface to reduce sample retention.


## Quality Control

Automated Quality Control procedures ensure that TipOne® is inspected for 100% performance. Every TipOne® product is certified RNase, DNase, DNA and Pyrogen free.

## Product Features

- **Graduated** at 2.5 µl and 10 µl
- **Universal fit:** One tip fits all! The tip collar is engineered to fit all popular pipettes.
- **Highly polished moulds:** For an extremely smooth surface to reduce sample retention.
- **Material:** All TipOne® tips are manufactured from pure, virgin polypropylene to ensure the highest quality product.
- **Certified:** Every TipOne® product is certified RNase, DNase, DNA and Pyrogen free.



10 µl Graduated TipOne® Tip, Racked (non-sterile)			
	Description	Pack Size	Cat. No.
	10 µl Graduated TipOne® Tip, Racked (non-sterile)	960 Tips (10 Racks × 96 Tips)	S1111-3800

Technical Specifications	
Autoclavable	Yes
Color Code Wafer	Red
Material	Polypropylene
Pack Size	960 Tips (10 Racks × 96 Tips)
Sterile	No
Color	Natural
Volume	10 µl
For use with	Abimed®, Biohit®, BRAND®, Capp®, Eppendorf®, Gilson®, Hamilton®, HTL, Nichiryo®, Rainin®, Sartorius®, Socorex®, STARLAB®, Thermo Scientific™