

100 µl UltraPoint® TipOne Graduated Filter Tip (sterile), Racked, Case



Nothing goes through!

Protect your samples against unwanted contamination by using TipOne® filter tips. The TipOne® filter is one of the best on the market and protects against aerosol cross-contamination without trapping your valuable samples. The filters contain no additives to interfere with samples.


During pipetting it is not possible to see with the eye if fine aerosols are going into your pipette. If they are, those aerosols can contaminate other samples at a later date.

TipOne® filter tips are proven to block >99% of aerosols and are among the best filters on the market. The results, shown right, speak for themselves!

**The comparative tests were conducted in an independent test institute with filter tips from other leading manufacturers.*

Product Features

- **Graduated** at 10 µl, 50 µl and 100 µl
- **Quality:** Highly effective HDPE filter
- **Safe:** Filter contains no additives to interfere with your samples
- **Retrieval:** The filter will not absorb samples, allowing recovery if necessary
- **Certified:** Sterile, and RNase, DNase, DNA and Pyrogen free
- **Sterile**

100 µl UltraPoint® TipOne Graduated Filter Tip (sterile), Racked, Case			
	Description	Pack Size	Cat. No.
	100 µl UltraPoint® TipOne Graduated Filter Tip (sterile), Racked, Case	7,680 Tips (80 Racks (case size) × 96 Tips)	S1123-1840-C

Technical Specifications	
Autoclavable	No
Color Code Wafer	Yellow
Material	Polypropylene
Pack Size	7,680 Tips (80 Racks (case size) × 96 Tips)
Sterile	Yes
Color	Natural
Volume	100 µl